



Healgen Covid-19 Antigen Test

The Coronavirus Ag Rapid Test Cassette (Swab) is an in vitro immunochromatographic assay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct nasopharyngeal (NP) swab specimens directly from individuals who are suspected of COVID-19 within the first ten days of symptom onset. This swab test is intended for the rapid detection of COVID-19 - results available within 15 minutes with high accuracy.

The Healgen Covid-19 Antigen Test is one of only three UK Government approved Rapid Antigen Test. It is in high demand from schools, councils and companies.

This test is being rolled out by the UK government to mass test the population in Liverpool and other cities, [click here](#) for UK Government approval (test 3 on the list).

- It is intended to aid in the rapid diagnosis of SARS-CoV-2 infections.
- The Healgen Covid-19 Antigen Test (Swab) does not differentiate between SARS-CoV and SARS-CoV-2.
- Each kit contains everything required to carry out the tests, no additional equipment is needed.
- Stage 3 approved by Public Health England & validated by the UK government.
- Results in 15 minutes.
- Each kit includes clear instructions and is easy to use.

Clinical Evaluation:

- Positive Percent Agreement (PPA): 97.3%
- Negative Percent Agreement (NPA): 99.6%
- Overall Percent Agreement (OPA): 98.9%
- Specimen: Nasopharyngeal (NP) Swab
- Time to Results: 15 minutes
- Shelf Life: usually 24 months from the date of manufacture


Benefits:

- Rapid testing for SARS-CoV-2 antigen within 15 minutes
- Facilitates patient treatment decisions quickly
- Simple, time-saving procedure
- All necessary reagents provided & no equipment needed
- High sensitivity and specificity

EACH KIT CONTAINS:

- 20 Test cassettes
- 20 Sterile swabs
- 20 Extraction tubes
- 1 Workstation
- 2 Buffers
- 1 Package insert & dropper tip

20 tests per box - 1000 tests per carton



Click here for a full instructional video on how to use the test kits